

pDRIVE-mROSA

A plasmid with the native mouse ROSA promoter

Catalog # pdrive-mrosa

For research use only

Version # 05H18-MT

PRODUCT INFORMATION

Content:

- 1 disk of lyophilized GT100 *E. coli* bacteria transformed by pDRIVE-mROSA GT100 genotype is: *F*-, *mcrA*, Δ (*mrr-hsdRMS-mcrBC*), Δ 80*lacZ* Δ M15, Δ *lacX74*, *recA1*, *endA1*.
- 4 pouches of *E. coli* Fast-Media® Zeo

Shipping and storage:

- Products are shipped at room temperature.
- Transformed bacteria should be stored at -20°C. Bacteria are stable up to one year when properly stored.
- Store *E. coli* Fast-Media® Zeo at room temperature. Fast-Media® is stable 18 months when stored properly.

Quality control:

- Plasmid construct has been confirmed by restriction analysis and sequencing.
- Bacteria have been lyophilized, and their viability upon resuspension has been verified.
- Promoter activity has been confirmed by transient transfection of 293 cells as well as other selected cell lines.

GENERAL PRODUCT USE

pDRIVE is an expression plasmid containing a native or composite promoter of interest. pDRIVE may be used to:

- **Subclone a promoter of interest into another vector.** Unique restriction sites are present at each end of the promoter allowing convenient excision. The 5' site is *Sda* I. *Sda* I is compatible with *Nsi* I and *Pst* I. The 3' restriction site is *Bsp*HI which includes the ATG start codon, and is compatible with *Nco* I, *Afl* III and *Sly* I.
- **Compare the activity of different promoters** in transient transfection experiments. Each pDRIVE promoter drives the expression of the *LacZ* reporter gene which allows for testing of the promoter's activity in transient transfection experiments. Furthermore, the *LacZ* gene is flanked by unique restriction sites (*Bsp*HI and *Eco*R I) for easy replacement with a different gene of interest.

PROMOTER CHARACTERISTICS

Element	Name	Origin	Size bp
Promoter	ROSA	Mouse	1926
5'UTR	ROSA	Mouse	457
Enhancer	-	-	-

The ROSA26 promoter, initially identified by random retroviral gene trapping in mouse embryonic stem cells¹, directs expression of reporter² and recombinase genes³ in all cells throughout embryonic development and in adult tissues. This TATA-less promoter is very effective *in vitro* in a very broad range of mammalian cell lines. The strength of the ROSA26 promoter is ascribed to the 10 potential Sp1 sites found within the CpG island extending from the proximal promoter to the first half of intron 1, the highest number of Sp1 sites ever recorded in any natural promoter. The mouse ROSA promoter provided by InvivoGen contains at its 3' end a synthetic intronic sequence featuring a consensus splice acceptor site.

1. Zambrowicz BP, Imamoto A. *et al.* 1997. Proc Natl Acad Sci USA. 94:3789-94
2. Kisseberth WC., Brettingen NT., Lohse JK., Sandgren EP. 1999. Dev Biol.214:128-38.
3. Farley FW, Soriano P, Steffen LS, Dymecki SM. 2000. Genesis. 28:106-10

PLASMID FEATURES

- **LacZ gene** encodes β -galactosidase an enzyme that catalyzes the hydrolysis of X-Gal, producing a blue precipitate that can be easily visualized under a microscope.
 - **SV40 pAn:** The Simian Virus 40 late polyadenylation signal enables efficient cleavage and polyadenylation reactions resulting in high levels of steady-state mRNA.
 - **pMB1 Ori** is a minimal *E. coli* origin of replication with the same activity as the longer Ori.
 - **EM7** is a bacterial promoter that enables the constitutive expression of the antibiotic resistance gene in *E. coli*.
 - **Sh ble** gene confers zeocin resistance therefore allowing the selection of transformed *E. coli* carrying a pDRIVE plasmid.
- Note: Stable transfection of clones cannot be performed due to the absence of an eukaryotic promoter upstream of the Sh ble gene.*

METHODS

Growth of pDRIVE-transformed bacteria:

Use sterile conditions to do the following:

- 1- Resuspend the lyophilized *E. coli* by adding 1 ml of LB medium in the tube containing the disk. Let sit for 5 minutes. Mix gently by inverting the tube several times.
- 2- Streak bacteria taken from this suspension on a zeocin LB agar plate prepared with the *E. coli* Fast-Media®Zeo agar provided (see below).
- 3- Place the plate in an incubator at 37°C overnight.
- 4- Isolate a single colony and grow the bacteria in TB supplemented with zeocin using the Fast-Media® Zeo liquid provided (see below).
- 5- Extract the pDRIVE plasmid DNA using the method of your choice.

Note: For long-term storage of the pDRIVE-transformed bacteria, prepare a 20% glycerol stock of the bacteria grown in the overnight liquid culture and freeze at -80°C.

Selection of bacteria with E. coli Fast-Media Zeo:

E. coli Fast-Media® Zeo is a **fast and convenient** way to prepare liquid and solid media for bacterial culture by using only a microwave. *E. coli* Fast-Media® Zeo is a TB (liquid) or LB (solid) based medium with zeocin. *E. coli* Fast-Media® Zeo can be ordered separately (reference # fas-zn-l, fas-zn-s).

Method:

- 1- Pour the contents of a pouch into a clean borosilicate glass bottle or flask.
- 2- Add 200 ml of distilled water to the flask
- 3- Heat in a microwave on MEDIUM power setting (about 400Watts), until bubbles start appearing (approximately 3 minutes). **Do not heat a closed container. Do not autoclave Fast-Media®.**
- 4- Swirl gently to mix the preparation. **Be careful, the bottle and media are hot, use heatproof pads or gloves and care when handling.**
- 5- Reheat the media for 30 seconds and gently swirl again. Repeat as necessary to completely dissolve the powder into solution. But be careful to avoid overboiling and volume loss.
- 6- Let agar medium cool to 45°C before pouring plates. Let liquid media cool to 37°C before seeding bacteria.


Note: Do not reheat solidified Fast-Media® as the antibiotic will be permanently destroyed by the procedure.

TECHNICAL SUPPORT

Toll free (US): 888-457-5873
Outside US: (+1) 858-457-5873
E-mail: info@invivogen.com
Website: www.invivogen.com

InvivoGen

3950 Sorrento Valley Blvd. Suite A
San Diego, CA 92121 - USA



175

SdaI (6)
1 CCTGCAGGTGAAGACGTTACACAAGTAACATGAGAAAGCAGAAAATGCAGGTCATCCACGCACCCCTGACCCAGGCCAGCAGGGCGGGCTGCAGCATCAG
101 TACACAGGAGAAAGATCCTTATTCTAAGAATGAGAAAGCAAAGGCGCCGATAGAATAAATTAGCATAGAAGGGGCTTTCCAGGAGTTAAAACCTTC

AvrII (251)
201 CTTCTGAGCGATTACCTACTAAAACCAGGGCTTTTGCCCACTACCATTACCTAGGATCTTGGCTTGACGGATTATAGGGGCATATCCCTCCCCCTCT

AvrII (337)
301 TCTTTAGAGTCGTTCTTAAAGATCGCTCTCCACGCCCTAGGCAGGAAAACGACAAAATCTGGCTCAATTCAGGCTAGAACCCTACAAATTCACAGG

EcoRV (403)
401 GATATCGCAAGGATACTGGGCATACGCCACAGGAGTCCAAGAATGTGAGGTGGGGGTGGCGAAGGTAATGTCTTTGGTGTGGAAAAAGCAGCAGCCAT
501 CTGAGATAGGAACTGAAAAACAGAGGAGAGCGCTTCAGGAAGATTATGGAGGGGAGGACTGGGCCCCACAGCGACCAGAGTTGTACAAGGCCGCAA
601 GAACAGGGGAGGTGGGGGCTCAGGGACAGAAAAAAGTATGTGTATTTGAGAGCAGGTTGGGAGGCCTCTCTGAAAAGGTATAAACGTGGAGTA

SalI (785)
701 GGCAATACCCAGGCAAAAAGGGGAGACCAGAGTAGGGGGAGGGGAAGAGTCTGACCCAGGGAAGACATTA AAAAGGTAGTGGGGTCTGACTAGATGAAGG
801 AGAGCCTTCTCTCTGGCAAGAGCGGTGCAATGGTGTGTAAGGTAGCTGAGAAGACGAAAAGGGCAAGCATCTTCTGCTACCAGGCTGGGGAGGCC
901 AGGCCACGACCCGAGGAGAGGGAACGAGGGAGACTGAGGTGACCCTTCTTCCCGGGGCCGGTCTGTGGTTCGGTGTCTTTTCTGTTGGAC

XhoI (1075)
1001 CTTACCTTGACCCAGGCGCTGCCGGGCTGGGCCGGGCTGCGGCGACGGCACTCCCGGAGGCAGCGAGACTCGAGTTAGGCCAACGCGCGGCCA
1101 CGCGCTTCTCTGGCCGGAATGCCCGTACCCTGAGGTGGGGTGGGGGCGAGAAAAGCGGAGCGAGCCCGAGGCGGGGAGGGGAGGCCAGGGGCC

NotI (1275)
1201 GAGGGGGCCGGCACTACTGTGTGGCGGACTGGCGGACTAGGGCTGCGTGAGTCTCTGAGCGAGCGGGCGCGGCCCTCCCGGGCGGGCGCA
1301 GCGGGCGGAGCGGGCGAGTCACTCAGCCCGTCCCGAGCGGAAACGCCACTGACCGCACGGGATTCCAGTGCCGGCGCCAGGGGCACGCGGGACA
1401 CGCCCCCTCCCGCGCGCCATTGGCTCTCCGCCACCGCCCCACACTTATTGGCCGGTGCGCCCAATCAGCGGAGGCTGCCGGGGCGCCTAAAGAA

NotI (1600)
1501 GAGGCTGTGCTTTGGGGCTCCGGCTCTCAGAGAGCCTCGGCTAGGtaggggatcgggactctggcgggagggcggttggtgctgttgcgggagtagggc
1601 ggccgcgagcagccctccgagcgtggtggagccgttctgtgagacagccgggtacgagtcgtgacgctggaaggggcaagcgggtggtgggagggaatgc

NcoI (1765)
1701 ggtccgacctgcagcaaccggaggggggagggagaagggagcggaaaagtctccaccggagcggccatggctcggggggggggggcagcggaggagcgc
1801 ttccggccagcgtctcgtcgtgattggcttctttctcccgctgtgtgaaaacacaattgtactaaccttcttctcttctctcctgacagGTGT

AgeI (1915) **BspHI (1931)** **NheI (1969)**
1901 GAAACAGGAAGAGAACCCGGTAGGAGGCCATCATGAGCGGTTCTCATCATCATCATCATGGTATGGCTAGCATGCTGGAGCAGCAAATGGGTGC
1 MetSer Gl ySer Hi sHi sHi sHi sHi sHi sHi sGl yMeAl aSer MeTThr Gl yGl yGl nGl nMeTGl yAr

Bsu36I (2030)
2001 GGATCTGTACGACGATGACGATAAGGTACCTAAGGATCAGCTTGAGGTTGATCCCGTCTTTTACAACGTCGTGACTGGGAAAACCTGGCGTTACCCAA
23 gAspLeuTyrAspAspAspLysValProLysAspGlnLeuGlyValAspProValValLeuGlnArgArgAspTrpGluAsnProGlyValThrGln
2101 CTTAATCGCCTTGACGACATCCCCCTTTCGCCAGCTGGCGTAATAGCGAAGAGGCCGACCAGTCCGCTTCCCAACAGTTGCGCAGCCTGAATGGCG
57 LeuAsnArgLeuAlaAlaHisProProPheAlaSerTrpArgAsnSerGluGluAlaArgThrAspArgProSerGlnLeuArgSerLeuAsnGlyG

Bsu36I (2267)
2201 AATGGCGCTTTGCTGGTTCCGGCACCAGAAGCGGTGCGGAAAGCTGGCTGGAGTGCATCTTCTGAGGCCGATACTGTGCTGCTCCCTCAAACCTG
90 IuTrpArgPheAlaTrpPheProAlaProGluAlaValProGluSerTrpLeuGluCysAspLeuProGluAlaAspThrValValValProSerAsnTr
2301 GCAGATGCACGGTTACGATGCGCCATCTACACCAACGTAACCTATCCATTACGGTCAATCCGCGTGTGTTCCACGGAGAAATCCGACGGGTTGTTAC
123 pGlnMetHisGlyTyrAspAlaProIleTyrThrAsnValThrTyrProIleThrValAsnProProPheValProThrGluAsnProThrGlyCysTyr
2401 TCGCTCATTAAATGTTGATGAAAGCTGGCTACAGGAAGGCCAGACGCAATTATTTTATGATGGCGTTAACTCGGCGTTTCATCTGTGGTGAACGGGC
157 SerLeuThrPheAsnValAspGluSerTrpLeuGlnGluGlyGlnThrArgIleIlePheAspGlyValAsnSerAlaPheHisLeuTrpCysAsnGlyAla
2501 GCTGGTTCGGTTACGGCCAGGACAGTCTGTTGCGCTCTGAATTTGACCTGAGCGCATTTTACGCGCGGAGAAAACCCGCTCGCGGTGATGGTGGC
190 rGTrpValGlyTyrGlyGlnAspSerArgLeuProSerGluPheAspLeuSerAlaPheLeuArgAlaGlyGlyAsnArgLeuAlaValMetValLeuAr
2601 TTGGAGTGACGGCAGTTATCTGGAAGATCAGGATATGTGGCGGATGAGCGGCATTTCCGTGACGCTCGTTGCTGCATAAACGACTACACAAATCAGC
223 gTrpSerAspGlySerTyrLeuGluAspGluAspMetTrpArgMetSerGlyIlePheArgAspValSerLeuLeuHisLysProThrThrGlnIleSer
2701 GATTTCCATGTTGCCACTCGCTTAATGATGATTTTCAGCGCGCTGTACTGGAGCTGAAGTTGAGTTGAGTTGAGTTGAGTTGAGTTGAGTTGAGTTGAGTT
257 AspPheHisValAlaThrArgPheAsnAspPheSerArgAlaValLeuGluAlaGlyValGlnMetCysGlyGlyLeuArgAspTyrLeuArgValT
2801 CAGTTTCTTTATGGCAGGTTGAAACGACGCTCGCCAGCGGCACCGCCCTTTCGGCGGTGAAATATCGATGAGCGTGGTGGTTATGCCGATCGGCTCAC
290 hrValSerLeuTrpGlnGlyGlyThrGlnValAlaSerGlyThrAlaProPheGlyGlyGlyIleIleAspGluArgGlyGlyTyrAlaAspArgValTh
2901 ACTACGCTGAACGTCGAAAACCCGAAACTGTGGAGCGCCGAAATCCGAATCTCTATCGTGGTGGTGGTGGTGGTGGTGGTGGTGGTGGTGGTGGTGGT
323 rLeuArgLeuAsnValGluAsnProLysLeuTrpSerAlaGlyIleProAsnLeuTyrArgAlaValValGluLeuHisThrAlaAspGlyThrLeuIle
3001 GAAGCAGAAGCTGCGATGTCGTTTCCGCGAGGTGCGGATTGAAAATGGTCTGCTGCTGCTGCTGCTGCTGCTGCTGCTGCTGCTGCTGCTGCTGCTGCT
357 GluAlaGluAlaCysAspValGlyPheArgGluValArgIleGluAsnGlyLeuLeuLeuLeuAsnGlyLysProLeuLeuIleArgGlyValAsnArgH

EcoRV (3156)
3101 ACGAGCATCATCTCTGCATGGTCAGGTCATGGATGAGCAGACGATGGTGCAGGATATCTGCTGATGAAGCAGAACTTTAACGCCGTGCGCTGTT
390 i sGluHi sHisProLeuHisGlyGlyValMetAspGluGlnThrMetValGlnAspIleLeuLeuMetLysGlnAsnAsnPheAsnAlaValArgCysSe

3201 GCATTATCCGAACATCCGCTGTGGTACACGCTGTGCGACCCTACGGCCTGTATGTGGTGGATGAAGCCAATATTGAAACCCACGGCATGGTGCCAATG
423▶ r Hi s Tyr P ro Asn Hi s P ro Leu T rp Tyr Thr Leu Cys Asp Arg Tyr Gl y Leu Tyr Val Val Asp Gl uAl a Asn I | e Gl u Thr Hi s Gl y Me t Val P ro Me t
3301 AATCGCTGACCGATGATCCGCCTGCCTACCGCGATGAGCGAACCGTAAACGCGAATGGTGCGACGGCATCGTAATACCCGAGTGATCATCTGGT
457▶ Asn Arg Leu Thr Asp Asp P ro Arg T rp Leu P ro Al a Me t Ser Gl u Arg Val Thr Arg Me t Val Gl n Arg Asp Arg Asn Hi s P ro Ser Val I | e I | e T rp S
3401 CGCTGGGAATGAATCAGGCACCGCGCTAATCAGCAGCAGCTGTATCGCTGGATCAAAATCTGTCGATCCTTCCCGCGGGTGAGTAAAGGCGGG
490▶ er Leu Gl y Asn Gl u Ser Gl y Hi s Gl y Al a Asn Hi s Asp Al a Leu Tyr Arg T rp I | e Lys Ser Val Asp P ro Ser Arg P ro Val Gl n Tyr Gl u Gl y Gl y
3501 AGCCGACACACGCGCACCGATATTATTTGCCGATGTACGCGCGCTGGATGAAGACCAGCCCTTCCCGCTGTGCCAAATGGTCCATCAAAAATGG
523▶ y Al a Asp Thr Thr Al a Thr Asp I | e I | e Cys P ro Me t Tyr Al a Arg Val Asp Gl u Asp Gl n P ro Phe P ro Al a Val P ro Lys T rp Ser I | e Lys Lys T rp
3601 CTTTCGCTACCTAGAGAGACGCGCCCGCTGATCCTTTGCCAATACGCCACGCGAATGGGTAACAGTCTGGCGTTCGCTAAATACCTGGCAGCGTTTC
557▶ Leu Ser Leu P ro Gl y Gl u Thr Arg P ro Leu I | e Leu Cys Gl u Tyr Al a Hi s Al a Me t Gl y Asn Ser Leu Gl y Gl y Phe Al a Lys Tyr T rp Gl n Al a Phe A
3701 GTCAATCCCCGTTACAGGGCGGCTTCGTCTGGACTGGGTGGATCAGTCCGTGATTAATATGATGAAACCGCAACCCGTTGGTCCGCTTACGGCGG
590▶ r g L n Tyr P ro Arg Leu Gl n Gl y Gl y Phe Val T rp Asp T rp Val Asp Gl n Ser Leu I | e Lys Tyr Asp Gl u Asn Gl y Asn P ro T rp Ser Al a Tyr Gl y Gl y
3801 TGATTTGGGATACCGCAACGCGCAATTCTGTATGAGCGCTGTCTGGTCTTGACCAGCGCCGCTTCCAGCCGCGGCGATTACGAGCAGCAACCGGAGCGG
623▶ y Asp Phe Gl y Asp Thr P ro Asn Asp Arg Gl n Phe Cys Me t Asn Gl y Leu Val Phe Al a Asp Arg Thr P ro Hi s P ro Al a Leu Thr Gl u Al a Lys Hi s Gl n
3901 CAGCAGTTTTCCAGTTCGGTTATCCGGGCAACCATCGAAGTGCACGCGAATACCTGTTCCGCTATAGCGATAACGAGCTCCTCGACTGGATGGTGG
657▶ Gl n Gl n Phe Phe Gl n Phe Arg Leu Ser Gl y Gl n Thr I | e Gl u Val Thr Ser Gl u Tyr Leu Phe Arg Hi s Ser Asp Asn Gl u Leu Leu Hi s T rp Me t Val A
4001 GGATTGATGGTAGTCGCTCCCGCGCGCTACCGCTGAGTGCCTGACCAAGGTAACAGTTGATTGAATTCAGAGCAGCGGAGCAAGCGGAGG
690▶ I a Leu Asp Gl y Lys P ro Leu Al a Ser Gl y Gl u Val P ro Leu Asp Val Al a P ro Gl n Gl y Lys Gl n Leu I | e Gl u Leu P ro Gl u Leu P ro Gl n P ro Gl u Se
4101 CGCCGGCAACTCTGGCTCACGTACCGTGTGCAACCAGCGCGATCGTCAAGCCGGGGCACATCAGCGCTGGCAGCAGTGGCGTCTGGCG
723▶ r Al a Gl y Gl n Leu T rp Leu Thr Val A rg Val Val Gl n P ro Asn Al a Thr Al a T rp Ser Gl u Al a Gl y Hi s I | e Ser Al a T rp Gl n Gl n T rp Arg Leu Al a
4201 GAAACCTCAGTGTGCGCTCCCGCGCTCCCGCATCCGCATCCGCTGGTGCCTGCACCAAGGTAACAGTTGATTGAATTTGATGAGCAGCGGAGCGGAGG
757▶ Gl u Asn Leu Ser Val Thr Leu P ro Al a Al a Ser Hi s Al a I | e P ro Hi s Leu Thr Thr Ser Gl u Me t Asp Phe Cys I | e Gl u Leu Gl y Asn Lys Arg T rp G
4301 AATTTAACCGCCAGTCCAGCTTCCTTTCACAGATGTGGATTGGCGATAAAAACAACTGCTGACGCCGCTGCAGTCACTCACCGTCCAGCTGGG
790▶ I n Phe Asn Arg Gl n Ser Gl y Phe Leu Ser Gl n Me t T rp I | e Gl y Asp Lys Lys Gl n Leu Leu Thr P ro Leu Arg Asp Gl n Phe Thr Arg Al a P ro Leu As
4401 TAACGACATTGGCGTAAGTGAAGCGACCGCATTTACCTAAGCGCTGGGTGAACGCTGGAAGCGCGGCGGCGATTACGAGCAGCGGAGCGGTTGTTG
823▶ p Asn Asp I | e Gl y Val Ser Gl u Al a Thr Arg I | e Asp P ro Asn Al a T rp Val Gl u Arg T rp Lys Al a Al a Gl y Hi s Tyr Gl n Al a Gl u Al a Al a Leu Leu
4501 CAGTGCACCGCATACACTTGTGTGCGGTGCTGATTACGACCGCTCACGCTGGCAGCATCAGGGGAAAACCTTATTATCAGCGGAAAACCTTACC
857▶ Gl n Cys Thr Al a Asp Thr Leu Al a Asp Al a Val Leu I | e Thr Thr Al a Hi s Al a T rp Gl n Hi s Gl n Gl y Lys Thr Leu Phe I | e Ser Arg Lys Thr Tyr A
4601 GGATTGATGGTAGTGTCAAATGGCGATTACCGTTGATTTGAAGTGGCGAGCGATACACCGCATCCGCGCGGATTGGCTGAACCTGCCAGCTGGCGCA
890▶ r gl | e Asp Gl y Ser Gl y Gl n Me t Al a I | e Thr Val Asp Val Gl u Val Al a Ser Asp Thr P ro Hi s P ro Al a Arg I | e Gl y Leu Asn Cys Gl n Leu Al a Gl
4701 GGTAGCAGAGCGGTAAACTGGCTCGGATAGGCGCGAAGAAAATATCCCGACCGCTTACTGCCCGCTGTTTTGACCGCTGGGTCTGCCATTGTCA
923▶ n Val Al a Gl u Arg Val Asn T rp Leu Gl y Leu Gl y P ro Gl n Gl u Asn Tyr P ro Asp Arg Leu Thr Al a Al a Cys Phe Asp Arg T rp Asp Leu P ro Leu Ser
BspLU11I (4802)
4801 GACATGTATACCCGTACGCTTTCCCGAGCGAAAACGGTCTGCGCTGCGGGACGCGGAATTGAATTATGGCCACACCAGTGGCGGGCGACTTCCAGT
957▶ Asp Me t Tyr Thr P ro Tyr Val Phe P ro Ser Gl u Asn Gl y Leu Arg Cys Gl y Thr Arg Gl u Leu Asn Tyr Gl y P ro Hi s Gl n T rp Arg Gl y Asp Phe Gl n P
4901 TCAACATCAGCCGTACAGTCAACAGCAACTGATGAAACGACCGTACCGCATCTGCTGACGCGGAAGAAGGCACATGGTGAATATCAGCGTTCCTCA
990▶ he Asn I | e Ser Arg Tyr Ser Gl n Gl n Gl n Leu Me t Gl u Thr Ser Hi s Arg Hi s Leu Leu Hi s Al a Gl u Gl u Gl y Thr T rp Leu Asn I | e Asp Gl y Phe Hi
5001 TATGGGGATTGGTGGCGGACTCCTGGAGCCGCTAGTATCGGCGAATACAGCTGAGCGCGGCTGCCTACCATTACAGTTGGTCTGGTGTCAAAA
1023▶ s Me t Gl y I | e Gl y Gl y Asp Asp Ser T rp Ser P ro Ser Val Ser Al a Gl u Leu Gl n Leu Ser Al a Gl y Arg Tyr Hi s Tyr Gl n Leu Val T rp Cys Gl n Lys
NheI (5122)
5101 TAATAATCTAGTCGAGAATTCGCTAGCTCGACATGATAAGATACATTGATGAGTTTGGACAAACCACAACAGTAACTAGAATGCAGTGAAAAAATGCTTTATTTG
1057▶ ●●●
5201 TGAATTTGTGATGCTATTGCTTTATTTGTGAAATTTGTGATGCTATTGCTTTATTTGTAACCATTATAAGCTGCAATAAACAAGTTAACAACAACAAATT
SwaI (5387)
5301 GCATTCATTTATGTTTCAGGTTCAAGGGGAGGTGTGGGAGTTTTTAAAGCAAGTAAACCTCTACAAATGTGGTAGATCCATTAAATGTTAATTA
5401 CTAGCCATGACCAAAATCCCTTAACTGAGTTTTCTGTTCCACTGAGCGTCAGACCCCGTAGAAAAGATCAAAGGATCTTCTTGTAGATCCTTTTTTCTGC
5501 GCGTAATCTGCTGTTGCAAAACAAAACCACCGCTACACGCGGTGGTGTGTTGCCGATCAAGAGCTACCAACTCTTTTTCCGAAGGTAACCTGGCTT
5601 CAGCAGAGCGCAGATACCAACTGTTCTTCTAGTGTAGCCGTAGTTAGGCCACCACTTCAAGAACTCTGTAGCACCCTACATACCTCGCTGTGTA
5701 ATCTGTTACCAAGTGGCTGCTGCCAGTGGCGATAAGTCTGTCTTACCGGGTTGGACTCAAGACGATAGTTACCGGATAAGGCGCAGCGGTGGGCTGAA
5801 CGGGGGTTCTGTGCACACAGCCAGCTTGGAGCGAACGACTACACCGAACTGAGATACCTACAGCGTGAGCTATGAGAAAAGCCACGCTTCCGAAGG
5901 GAGAAAAGCGGACAGGTATCCGTAAGCGGCGAGGTCGGAACAGGAGAGCGACGAGGGAGCTTCCAGGGGAAACGCCTGGTATCTTTATAGCTCTGTC
6001 GGGTTTCGCCACCTCTGACTTAGCGTCTGATTTTTGTGATGCTCTGAGGGGCGGAGCCTATGAAAAACGCCAGCAACCGCCCTTTTTACGGTTC
BspLU11I (6125) AseI (6163)
6101 TGGCCTTTTGTGCTGCTGCTTTGTCTCACATGTTCTTAATAAATTTTCAAAGTAGTTGACAATTAATCATCGGCATAGTATATCGGCATAGTATAATACG
6201 ACTCACTATAAGGAGGCCATCATGGCCAAGTTGACCACTGCTGTCCAGTGTCTCAGCGCAGGATGTGGCTGGAGCTGTTGAGTTCTGGACTGACAGGT
1▶ Me t Al a Lys Leu Thr Ser Al a Val P ro Val I Leu Thr Al a Arg Asp Val Al a Gl y Al a Val Gl u Phe T rp Thr Asp Arg L
6301 TGGGGTTCTCCAGAGATTTGTGGAGGATGACTTTGCAAGGTGTTGCTGAGAGATGATCACCCTGTCTCAGCAGTCCAGGACAGGTTGGCTGA
27▶ eu Gl y Phe Thr Arg Asp Phe Val Gl u Asp Asp Phe Al a Gl y Val Al a Arg Asp Val I Leu Thr Phe I | e Ser Al a Val Gl n Asp Gl n Val I P ro As
6401 CAACACCTGGCTTGGGTGTGGGTGAGAGGACTGGATGAGCTGTATGCTGAGTGGAGTGGTGGTCTCCCAACTTCCAGGATGCCAGTGGCCCTGCC
60▶ p Asn Thr Leu Al a T rp Val T rp Val A rg Gl y Leu Asp Gl u Leu Tyr Al a Gl u T rp Ser Gl u Val Val Ser Thr Asn Phe Arg Asp Al a Ser Gl y P ro Al a
6501 ATGACAGAGATTGGAGAGCAGCCCTGGGGGAGAGTHTGCCCTGAGAGACCCAGCAGCACTGTGTGCACTTTTGTGGCAGAGGAGCAGGACTGAGGAT
94▶ Me t Thr Gl u I | e Gl y Gl n P ro T rp Gl y A rg Gl u Phe Al a Leu Arg Asp P ro Al a Gl y Asn Cys Val Hi s Phe Val Al a Gl u Gl n Asp ●●●
6601 AAGAATGAGTTTCAGAAAAGGGGCGCTGAGTGGCCCTTTTTCAACTAATTA